

Practitioner Testimonials

As a qualified practitioner, I am always looking for the best products to recommend to my clients. This new ionic formulation silver is exciting.

I believe this product to be the best value, cleanest and most effective on the market. I will urge my clients to insist on using this 21st century technology and throw away the colloidal silvers of yesteryear.

Ms. Dawn Durbin N.D.,
Dip. Remedial Therapy, Dip. Of Counselling

Part of my professional approach as a Naturopath is to identify clients' exposure to heavy metals and toxins, which could be affecting their health.

I am also very proactive in dispensing supplements that are ionic in nature and form. Hence, when HealthWest Ionic silver came onto the market, I have been able to recommend it with confidence. It is clear, tasteless, and odourless. I have also educated clients on the use of ionic silver for their animals – both farm and domestic. I am impressed with this company's high standards of manufacturing strategies to ensure preservation of the environment. Their approach is truly holistic.

Mavis Carruthers ND

Disclaimer

No therapeutic claims are made and no medical advice is offered. If you have a health problem we recommend seeking medical advice from your qualified health professional. This information is strictly a source of general information and is not intended for use as a tool for self-diagnosis.

How to take Ionic Silver

Ionic Silver can be consumed undiluted or mixed with good quality water.

Can be taken orally, applied topically, gargled or nebulised.

When taking Ionic Silver orally, drink 3-4 glasses of good quality water daily to help flush out any infective agents.

Acute conditions

Adults 20-30ml hourly

Children under 12 10-15ml hourly

When condition improves take for a further three days at:

Adults 20 to 30ml twice daily

Children under 12 10ml twice daily

As a preventative

Adults 10-20ml twice daily

Children under 12 10ml twice daily

Dogs and Cats

Hourly in acute conditions

Large dogs, over 45kg 30ml

Dogs 20 to 40kg 20ml

Dogs 5 to 20kg 15ml

Small dogs and cats up to 5kg 10ml

Puppies and kittens 5ml

Also Available



Health Products That Work

Your Local Stockist

PURE

NO METAL

IONIC SILVER

Silver Ions have the unique ability to disable micro-organisms and render them harmless



www.healthwest.com.au

Product of Australia

The Truth about Colloidal Silver & Ionic Silver

All colloidal silver contains particles of metal suspended in distilled water

HealthWest Ionic Silver does not contain any metal


HealthWest Pure Ionic Silver is deliberately designed to entirely overcome the inherent limitations of colloidal silver and provide a pure Silver Ion solution completely free of metal.

It has a built-in delivery mechanism to shield it from chloride and other salts in the body, enabling it to be delivered intact to where it is needed.

HealthWest Ionic Silver... How does it work?

Microbes, like all living organisms rely on an electrical or nervous system to sustain life, and this is how they live communicate and breed.


How can you tell if there is metal in a solution?



Charged silver ions disable this internal electrical system and the microbe ceases to function and simply collapses and dies.

Because of the unique way that silver ions work there can be no resistance either inherited or developed.

HealthWest complies with the Australia Food Safety Standard and our products are analyzed to ensure their purity from a NATA certified testing laboratory.



“ Research by Rice University in the USA found that it is the silver ions, not the particles themselves that are toxic to bacteria. ”

In other words, it is the charged silver ions that kill germs not the metal particles or colloids.

TRUTH

- All Colloidal Silver contains metal
- HW Ionic Silver does not contain metal
- Silver ions have a positive electrical charge
- Colloidal silver contains colloids which are electrically neutral.
- It is the silver ions not the colloids that kill germs

MYTH

- The smaller the particle size, the more effective.
- Particles or colloids have positive electrical charge.

Normally when a laser is shone through a liquid other than distilled water, a path of light can be seen similar to when you shine a torch through smoke or fog.

All colloidal silvers will show this path of light as the laser beam highlights the metal particles contained within the liquid.

This test will clearly demonstrate that HealthWest Pure Ionic Silver contains only Silver Ions and does not contain any colloids, metals or nanoparticles.

These Silver Ions are made with a built-in delivery mechanism to shield them from chloride and other salts, enabling them to be delivered intact to where they are needed.

All this is achieved without any metal or particles in the solution as is evidenced by the Tyndall test.

NEW
PRODUCTS

Pure Ionic Silver Practitioner Strength

HealthWest Pure Ionic Silver is deliberately designed to entirely overcome the inherent limitations of colloidal silver and provide a pure Silver Ion solution completely free of metal.

It has a built-in delivery mechanism to shield it from chloride and other salts in the body, enabling it to be delivered intact to where it is needed.

**HealthWest
Ionic Silver...
How does it work?**

Microbes, like all living organisms rely on an electrical or nervous system to sustain life, and this is how they live communicate and breed.

Charged silver ions disable this internal electrical system and the microbe ceases to function and simply collapses and dies.

Because of the unique way that silver ions work there can be no resistance either inherited or developed.

HealthWest complies with the Australia Food Safety Standard and our products are analyzed to ensure their purity from a NATA certified testing laboratory.



NEW
PRODUCTS

Pure Ionic Silver Spray and Spray Refill Spray for instant relief Does not sting

Silver ions are increasingly being recognized for their broad-spectrum antimicrobial qualities and bandages impregnated with silver ions are especially in demand. They prevent bacterial growth and therefore speed healing time, making them extremely valuable for treating burn and wound victims. The fact that silver ions are effective against a very broad range of bacteria is well established and due to recent advances in the delivery of silver ions, they are being used in a rapidly growing range of wound care, medical and industrial products.



Ear Nose & Throat Spray Extra Strength ENT PS50

Spray as often as desired

HealthWest Ear Nose & Throat spray is designed to get at those hard to reach places, and delivers a spray of positively charged silver ions and pure water.

Does not contain any nanoparticles, (metal colloids), fillers, stabilisers or chemicals of any sort.